# 1.10003.0001



# Cu

# 1. Method

Copper(II) ions are reduced to copper(I) ions by a reducing-agent mixture. The copper(I) ions react with 2,2'-biquinoline (cuproin) to form a violet complex. The copper concentration is measured **semi-quantitatively** by visual comparison of the reaction zone of the test strip with the fields of a color scale.

# 2. Measuring range and number of determinations

Measuring range / color-	Number of	
scale graduation	determinations	
<b>10</b> - 30 - 100 - <b>300 mg/l Cu</b>	100	

# 3. Applications

This test measures both copper(II) and copper(I) ions. It is also suited for the detection of copper in metallic materials and surfaces (see section 7).

#### Sample material:

Drinking water

Swimming-pool water

Wastewater, especially from the electroplating and printing industries Allovs

Wine

# 4. Influence of foreign substances

This was checked in solutions with 30 and 0 mg/l Cu. The determination is not yet interfered with up to the concentrations of foreign substances given in the table.

Concentrations of foreign substances in mg/l					
Ag+	1000	Fe <sup>2+</sup>	1000	NH₄+	1000
Al <sup>3+</sup>	1000	Fe <sup>3+</sup>	1000	Ni <sup>2+</sup>	1000
Ba <sup>2+</sup>	1000	[Fe(CN) <sub>6</sub> ] <sup>4-</sup>	1	NO <sub>2</sub> -	1000
Ca <sup>2+</sup>	1000	[Fe(CN) <sub>6</sub> ] <sup>3-</sup>	1	NO <sub>3</sub> -	1000
Cd <sup>2+</sup>	1000	F .	250	Pb <sup>2+</sup>	1000
Cl	1000	K⁺	1000	PO₄ <sup>3-</sup>	1000
CN-	1	Mg <sup>2+</sup>	1000	SO32-	1000
Co <sup>2+</sup>	1000	MnO <sub>4</sub> -	1000	SO42-	1000
CrO42-	500	Na⁺	1000	Zn <sup>2+</sup>	1000

# 5. Reagents and auxiliaries

The test strips are stable up to the date stated on the pack when stored closed at +15 to +25  $^\circ$ C.

#### Package contents:

Tube containing 100 test strips

#### Other reagents:

MColorpHast<sup>™</sup> Universal indicator strips pH 0 - 14, Cat. No. 109535

Sodium acetate anhydrous for analysis EMSURE®, Cat. No. 106268

Sulfuric acid 0.5 mol/l TitriPUR®, Cat. No. 109072 Copper standard solution CertiPUR®, 1000 mg/l Cu, Cat. No. 119786

## 6. Preparation

- Samples containing more than 300 mg/l Cu must be diluted with distilled water.
- The pH must be within the range 2 6. If necessary, buffer the sample with sodium acetate or, respectively, adjust the pH with sulfuric acid.

# 7. Procedure

Immerse the reaction zone of the test strip in the pretreated sample (15 - 25 °C) for 1 sec. Shake off excess liquid from the strip and after 30 sec

Shake off excess liquid from the strip and **after 30 sec** determine with which color field on the label the color of the reaction zone coincides most exactly. Read off the corresponding result in mg/l Cu.

#### Notes on the measurement:

- The color of the reaction zone may continue to change after the specified reaction time has elapsed. This must not be considered in the measurement.
- If the color of the reaction zone is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 300 mg/l Cu is obtained.

Concerning the result of the analysis, the dilution (see also section 6) must be taken into account:

Result of analysis = measurement value x dilution factor

#### Determination on metallic surfaces:

Moisten the reaction zone of the test strip with water and lightly press it on the surface to be tested (e.g. that of a coin) for approx. 10 - 30 sec. The obtained measurement results are only guideline values. A violet color is produced by an amount of copper as low as 0.5  $\mu$ g.

# 8. Method control

To check test strips and handling: Dilute the copper standard solution with distilled water to 100 mg/l Cu and analyze as described in section 7.

Additional notes see under www.qa-test-kits.com.

## 9. Note

Reclose the tube containing the test strips immediately after use.



EMD Millipore Corporation, 290 Concord Road, Billerica, MA 01821, USA, Tel. +1-978-715-1335