

1.14598.0001
1.14598.0002Spectroquant®
Fluoride TestF⁻

1. Method

In a buffered, weakly acidic solution, fluoride ions react with alizarin complexone and lanthanum(III) to form a violet complex that is determined photometrically. The method is analogous to EPA 340.3 and APHA 4500-F E.

2. Measuring range and number of determinations

Cell mm	Measuring range mg/l F ⁻	Number of determinations
10	0.10 - 2.00 1.0 - 20.0	100 (Cat. No. 1.14598.0001) or 250 (Cat. No. 1.14598.0002)

For programming data for selected photometers / spectrophotometers see www.service-test-kits.com.

3. Applications

This test does not measure any complex fluorine compounds as occur in some cases in galvanic waste water.

Sample material:

Groundwater and surface water, sea water
Drinking water and mineral water
Wastewater and percolating water

4. Influence of foreign substances

This was checked in solutions containing 1 (10) and 0 mg/l F⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. The values in parentheses apply for the measuring range 1.0 - 20.0 mg/l F⁻.

Concentrations of foreign substances in mg/l or %			
Ag⁺	1 (10)	Mn²⁺	2 (20)
Al³⁺	1 (5)	NH₄⁺	1000 (1000)
Ca²⁺	50 (500)	Ni²⁺	1 (10)
Cd²⁺	5 (50)	NO₂⁻	1000 (1000)
CN⁻	0.5 (5)	Pb²⁺	1 (10)
CO₃²⁻	1000 (1000)	PO₄³⁻	2 (20)
Cr³⁺	5 (50)	S²⁻	100 (1000)
Cr₂O₇²⁻	1000 (1000)	SiO₃²⁻	1000 (1000)
Cu²⁺	2 (25)	Zn²⁺	2 (25)
Fe³⁺	2 (20)	EDTA	0.1 (1)
Hg²⁺	2 (20)	Free chlorine	10 (50)
Mg²⁺	100 (1000)	Anionic surfactants ¹⁾	1000 (1000)
		Cationic surfactants ²⁾	50 (500)
		Nonionic surfactants ³⁾	1000 (1000)
		Na-acetate	10 % (20 %)
		NaCl	10 % (20 %)
		NaNO ₃	10 % (20 %)
		Na ₂ SO ₄	10 % (20 %)

¹⁾ tested with Na-dodecyl sulfate

²⁾ tested with N-cetyl-N,N,N-trimethylammonium bromide

³⁾ tested with Triton[®] X-100

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

Reagent F-1: 1 bottle

Reagent F-2: 1 bottle (Cat. No. 1.14598.0001) or
2 bottles (Cat. No. 1.14598.0002)

2 AutoSelectors

Other reagents and accessories:

MColorpHast™ Universal indicator strips pH 0 - 14, Cat. No. 109535

Sodium hydroxide solution 1 mol/l TitriPUR®, Cat. No. 109137

Sulfuric acid 0.5 mol/l TitriPUR®, Cat. No. 109072

Water for analysis EMSURE®, Cat. No. 116754

Fluoride standard solution CertiPUR®, 1000 mg/l F⁻, Cat. No. 119814

Pipettes for pipetting volumes of 0.50 and 5.0 ml

Pipette for a pipetting volume of 2.0 ml - suited for pipetting volatile organic solvents

Rectangular cells 10 mm (2 pcs), Cat. No. 114946

6. Preparation

● The pH must be within the range 3 - 8.

Adjust, if necessary, with sodium hydroxide solution or sulfuric acid.

● Filter turbid samples.

7. Procedure

7.1 Measuring range 0.10 - 2.00 mg/l F⁻

Reagent F-1	2.0 ml	Pipette into a test tube.
Pretreated sample (10 - 40 °C)	5.0 ml	Add with pipette and mix.
Reagent F-2	1 level blue microspoon (in the cap of the F-2 bottle)	Add and shake vigorously until the reagent is completely dissolved.
Leave to stand for 5 min (reaction time) , then fill the sample into a 10-mm cell, and measure in the photometer.		

7.2 Measuring range 1.0 - 20.0 mg/l F⁻

Reagent F-1	2.0 ml	Pipette into a test tube.
Distilled water ¹⁾ (10 - 40 °C)	5.0 ml	Add with pipette and mix.
Pretreated sample (10 - 40 °C)	0.50 ml	Add with pipette and mix.
Reagent F-2	1 level blue microspoon (in the cap of the F-2 bottle)	Add and shake vigorously until the reagent is completely dissolved.
Leave to stand for 5 min (reaction time) , then fill the sample into a 10-mm cell, and measure in the photometer.		

¹⁾ It is recommended to use water for analysis EMSURE®, Cat. No. 116754.

Notes on the measurement:

- **Reclose the reagent bottles immediately after use.**
- **Certain photometers may require a blank** (preparation as per measurement sample, but with distilled water instead of sample).
- For photometric measurement the cells must be clean. Wipe, if necessary, with a clean dry cloth.
- Measurement of turbid solutions yields false-high readings.
- The pH of the measurement solution must be within the range 4.4 - 4.7.
- The color of the measurement solution remains stable for 60 min after the end of the reaction time stated above.
- In the event of fluoride concentrations exceeding 5 mg/l (for the procedure described in section 7.1) or exceeding 50 mg/l (for the procedure described in section 7.2), other reaction products are formed and false-low readings are yielded. In such cases it is advisable to conduct a plausibility check of the measurement results by diluting the sample (1:2, 1:10).

8. Analytical quality assurance

recommended before each measurement series

To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working, a dilute fluoride standard solution containing 1.00 or 10.0 mg/l F⁻ can be used.

Sample-dependent interferences (matrix effects) can be determined by means of standard addition.

Additional notes see under www.qa-test-kits.com.

Characteristic quality data:

In the production control, the following data were determined in accordance with ISO 8466-1 and DIN 38402 A51:

	Measuring range mg/l F ⁻	
	0.10 - 2.00	1.0 - 20.0
Standard deviation of the method (mg/l F ⁻)	± 0.015	± 0.15
Coefficient of variation of the method (%)	± 1.4	± 1.4
Confidence interval (mg/l F ⁻)	± 0.04	± 0.4
Number of lots	30	30

Characteristic data of the procedure:

	Measuring range mg/l F ⁻	
	0.10 - 2.00	1.0 - 20.0
Sensitivity: Absorbance 0.010 A corresponds to (mg/l F ⁻)	0.02	0.2
Accuracy of a measurement value (mg/l F ⁻)	max. ± 0.12	max. ± 1.2

For quality and batch certificates for Spectroquant® test kits see the website.

9. Note

Information on disposal can be obtained at www.disposal-test-kits.com.

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.analytical-test-kits.com

EMD Millipore Corporation, 290 Concord Road,
Billerica, MA 01821, USA, Tel. +1-978-715-4321

