SC-110 Conductivity/ Resistivity Portable Meter

Operation Manual



CONTENTS

1.	Specifica	tions	1
2.	Configur	ation	
	2.1 Displa	ay	2
	2.2 Socke		2
	2.3 Keypa	ad	3
3.	Operatio	n	
	3.1 Meas	urement	4
	3.2 Meas	urement parameter and measurement range	4
	3.3 Temp	erature function and compensation	5
	3.4 Auto	read	5
4.	Settings		
	Block dia	gram of settings	6
	4.1 Parar	neter settings	7
	4.1.1	Entry of set-up mode	7
	4.1.2	Settings of measurement mode	7
	4.1.3	Digital filter settings	7
	4.1.4	Data overwriting function	8
	4.1.5	Settings of RS-232 digital communication	9
	4.1.6	System time settings	10
	4.1.7	Auto power-off settings	11
	4.2 Delete	e all of the data	
	4.3 System	m Reset	
	Plaak dia	grow of solibustion and temperature companyation settings	12
	A A Sottin	gr of colibration parameter	13 14
	4.4 Settin	Tomporature componentian coefficient setting for standard solv	¹⁴
	4.4.1	calibration	14
	442	Multi-noint of standard solution calibration setting	14
	443	Temperature compensation coefficient setting for measurement	t mode 15
	т.т.5	remperature compensation coefficient setting for measurement	1110uc_15
	Block dia	gram of data storage and transmitting	16
	4.5 Data	storage	17
	4.5.1	Time-interval auto data storage	17
	4.5.2	Manual data storage	18
	4.5.3	Overflow of data storage	18

4.6 Data transmitting and display		
4.6.1	Time-interval auto transmitter	19
4.6.2	Data display & manual transmitting	20

5.	Calibrati	on	
	Block dia	ngram of calibration	25
	5.1 Stand	ard solution preparation	26
	5.2 Calib	ration mode	26
	5.3 Entry	of calibration mode	26
	5.3.1	Cell constant calibration	26
	5.3.2	Standard solution single-point calibration	26
	5.3.3	Standard solution five-point calibration	27
	5.3.4	Error message	27
6.	Introduc	tion set for RS-232	28
7.	Modbus	protocol and instructions for SC-110	29
8.	Optional	accessories	38
9.	Trouble s	shooting	35
10.	Maintena	ance	40

1. Specifications

Mode		SC-110		
Measuring Pa	arameter	Conductivity / Resistivity / Salinity / TDS / Temp.		
	Cond.	0.00µS/cm ~ 200.0 mS/cm 5 measuring ranges for manual or auto range		
	Res.	0.00~100.00 MΩ·cm		
Range	Salinity	0.0~70.0 ppt		
	TDS	0~2000 ppm		
	Temp.	-10.0~110.0°C		
	Cond.	0.01 µS/cm		
	Res.	0.01 MΩ·cm		
Resolution	Salinity	0.1 ppt		
	TDS	1 ppm		
	Temp.	0.1°C		
	Cond.	±0.5%(±1Digit)		
	Res.	$\pm 1\%(\pm 1$ Digit)		
Accuracy	Salt.	±0.5%(±1Digit)		
	TDS	±0.5%(±1Digit)		
	Temp.	$\pm 0.2^{\circ}C(\pm 1Digit)$		
Temperature Compensation		PT-1000 or NTC 30K auto recognition		
		Auto / Manual selectable		
Tempera Coeffici	nture ient	Non-linear / Linear (compensation from $0.00 \sim 10.00$ %)		
Ref. Tempe	erature	T _{ref} 25°C		
Manual cali	bration	Single point known standard solution or cell constant calibration		
Standard s calibrat	olution tion	Single or up to 5-point known standard solution calibration; (10 / 84 / 1,413 µS/cm, 12.88 / 100 mS/cm)		
RS232 int	erface	Yes		
Output Data	logging	450 sets (measuring value, temp., date, time & ID)		
Ambient temperature		0~50°C		
Storage Temperature		-20~70°C		
Display		LCD Display		
Power		4 x AA battery		
Battery		Around 350 ~ 1000 hours according to measuring range Programmable auto shut-down		
Protect	ion	IP65		
Dimens	ion	185mm x 98mm x 38mm (L x W x H)		
Weigl	ht	0.35Kg		

2. Configuration

2.1 Display:



2.3 Keypad:



Mode

-

Cal

Auto read

Ūp

Down

Enter

Recall

- : Power switch
- : Conductivity/Salinity/TDS measurement shift key; press it at any time will allow back to the measurement mode
- : Enter calibration mode
- : Key for auto read in measurement mode. Press one time for startup and press again for cancellation
- : Increase value / Scroll upward or leftward
- : Decrease value / Scroll downward or rightward
- : Confirm input / Execute / Transmitting data manually via RS-232 in measuring mode
- : Manual data display, and roll up for next data transmitting function (See Chapter 4.6)



: Manual data storage



: Enter system parameter settings. (<u>Cond.</u> / <u>Res.</u> mode selection, auto turn-off, etc) (See Chapter 4.1)



: Delete all of stored data (See Chapter 4.2)



: System Reset (See Chapter 4.3)



: Shift among the five conductivity measurement manual range & auto range (See Chapter 3.2)



: Enter setting of calibration (See Chapter 4.4)



: Enter time-interval auto data storage menu (See Chapter 4.5.1)



: Enter time-interval auto data transmitting (See Chapter 4.6.1)

3. Operation

3.1 Measurement:

Press **POWER** key to start up the instrument. It will automatically enter the measurement mode of final operation and begin to measure.

3.2 Measurement unit and measurement range switch

- 1. Press **MODE** key to switch among Conductivity, Salinity, and TDS in <u>Cond.</u> Mode. This function is unavailable in <u>Res.</u> mode.
- 2. Press **ENTER** + **MODE** keys to shift among the five manual ranges and the auto-range. After entering auto mode, AUTO characters will disappear in 3 seconds.



Remark: Before measuring salt and TDS, please finish the calibration according to the conductivity calibration steps, and then set the factor value according to the standard solution of salt and TDS.

3.3 Temperature function and compensation

- ATC: Auto judgment of whether there is temperature probe and recognition of NTC30K or PT1000. Press UP or DOWN key to adjust the temperature within ±5.0°C.
 Press UP + DOWN key simultaneously to back to the original value.
- MTC: Enter MTC 25°C automatically when there is no temperature probe. Press UP key or DOWN key to adjust the temperature. Press UP + DOWN key simultaneously to back to the default value 25.0°C.

3. **Temperature compensation:**

Because the temperature of the sample solution may not be 25°C, for standardization, the conductivity value of solution not at the temperature of 25°C will be calculated to that of solution at the temperature of 25°C. This is temperature compensation.

- 4. According to the temperature coefficient required by the measurement, three kinds of temperature compensation can be applied: Lin, nLin, and non compensation (0.00%). See the setting method in 4.3.3, and the user can set as needed.
- 5. **Temp. Coefficient** (called as TC for short), the conductivity value of the solution will increase as the temperature rises at a different rate. See the relations as follows:

C25	Conductivity at 25°C	$Formula 1: Ct = Cor (1+\beta(T-25))$
Ct	Conductivity at T°C	Formula 1. Ct = C25 { $1 \pm p(1 - 23)$ }
Т	Temp. of the sample solution	Example 2: $\theta = (Ct - C_{22}) / (C_{22} (T - 25))$
β	Temp. Coefficient	Formula 2. $p = (C1 - C25) / \{C25 (1 - 23)\}$

6. How to measure the TC of the tested solution: If the user wants to get a more accurate value, the above formula can be applied to measure the TC of the tested solution and set an appropriate TC on the machine. Take 0.01N KCl as an instance, set the TC as non compensation (0.00%) in the machine. Control the same tested solution respectively at 5°C and 20°C. The value 1413 μ S tested at C25 is just the conductivity of the solution at 25°C; and the value 1278 μ S tested at Ct is just the conductivity of the solution at 20°C. According to formula 2, $\beta=1.91\%$.

3.4 Auto Read:

Under measurement mode, press **AUTO READ** key to hold the present reading(HOLD), and press **ENTER** key to confirm the Auto-Read function and to activate new measurement. Wait until the measurement is stable, the instrument holds the value automatically(HOLD) to finish the procedure. To cancel the hold value, press **AUTO READ** key or **MODE** key to relieve the Auto-Read hold value.



4. Settings

Block diagram of settings



4.1 Parameter settings:

4.1.1 Entry of set-up mode:

While the instrument's power is off, press **MODE** key and hold the pressing continuously, and press **POWER** key until the display shows model of the instrument. Release the **POWER** key and wait until the menu display enter into parameter settings mode, and then release **MODE** key. Press upward or downward key to select setting menu. To exit the parameter settings mode, press **MODE** key to back to measurement mode.

4.1.2 Settings of measurement mode:



4.1.3 Digital filter settings (Sample average of measurements)

Enter the setup of Digital filter. You may select the number of sample to be averaged each time to become a reading which is gradually counted in order to increase the stability of measurement.



4.1.4 Data overwriting function

Enter "Data overwriting function" settings



4.1.5 Settings of RS-232 digital communication:

Enter "Settings of RS-232"



4.1.6 System time settings

Users may set the year, month, date, hour, and minute. When replacing the batteries, it is necessary to enter to the setting to update the system time.



4.1.7 Auto power-off settings

The auto power off timer setting range: 10 minutes minimum, 23 hours and 59 minutes maximum.



4.2 Delete all of the data

Under power off status, press & hold **STORE** key and then press & hold **POWER** key until the display shows the model of instrument. Release **POWER** key. While the menu display shows data clean text, release **STORE** key.

4.3 System Reset

In the power-off status, push **CAL** first and then **POWER**; after that, release **POWER** first and then release **CAL** to enter the system reset page.



Factory defaults:

Parameter: Conductivity Digital filter: 6 Data over-writing: Auto RS-232: ASCII, non, 9600 System time: 2000, 1/1, 00:00 Auto power-off: Auto, 10 minutes

Block diagram of calibration and temperature compensation settings



4.4 Settings of calibration parameter:

In the power-on status, press and hold **ENTER** key first and then press **CAL** key to enter set-up mode of calibration parameter. Use upward key or downward key to select the menu.

4.4.1 Temperature compensation coefficient setting for cell constant calibration

(Only suitable when the measurement unit is conductivity)



4.4.2 Multi-point of standard solution calibration setting

(Only suitable when the measurement unit is conductivity)



4.4.3 Temperature compensation coefficient setting for measurement mode





Block diagram of data storage and transmitting

4.5 Data storage

Under data storage mode, users may store conductivity or resistivity or TDS or Salinity measurement with both No.(numbering) and ID (two kinds of classification). The "Number." represents the order of storage sequence. When data is stored to the 450^{th} data, the display will shows "Full" to represent there is no more space for data logger memory. While ID category $1\sim1000$ is for grouping the samples which the user may require. The Auto data storage can be set according to the time interval setting, the minimum interval is 5 seconds. If turn on the Data Rewritable Function, when the data memory is full the data logger will overwrite the new data.

4.5.1 Time-interval auto data storage

Under the measurement mode, press and hold **ENTER** key and press **STORE** key simultaneously to enter the auto storage setting.



4.5.2 Manual data storage

In the measurement mode, press **STORE** key to enter manual store page. Press **MODE** key to exit at any time.



4.5.3 Overflow of data storage

When the data is full to store, the display will shine as a reminder, and five seconds later, it will back to the original measurement mode.

Users may delete all of the data according to the ch 4.2 instruction by pressing **STORE** + **POWER** key while power off. (When auto rewritable function is turned on, there is no this data full display reminder.)

4.6 Data transmitting and display

4.6.1 Time interval auto transmitting

Transmit the real-time measurement to a computer or a printer through RS-232) Under measurement mode, press and hold **ENTER** key and then press **RECALL** key to enter auto transmitting menu display. After finishing setting, the data is transmitted through RS-232 interface.



4.6.2 Data display & manual transmitting

In the measurement mode, press **RECALL** key to roll up the data display & manual transmitting menu, and press **ENTER** to enter.

<u>1. Display the stored data</u>



2. Manual stored data transmitting

Transmit all stored data through RS-232 In the measurement mode, press **RECALL** 2 times to find the menu.



3. Display the last calibration data

In the measurement mode, press **RECALL** key 3 times to find the menu.

(1). Read the CELL constant calibration data

(when the calibration parameter is set as CELL constant calibration)



(2). Read single-point calibration data of standard solution (when the calibration parameter is set as Std 1P calibration)



(3). Read five-point calibration data of standard solution (when the calibration parameter is set as Std 5P calibration).



4. Calibration data transmitting

Transmit the stored cal. data through RS-232. In the measurement mode, press **RECALL** 4 times to find the menu.



5. Calibration

Block diagram of Calibration



5.1 Standard solution preparation:

Put the reagent potassium chloride in an oven of 150° C ~ 180° C for 5 hours, and then put it in a dry container until it cools to the room temperature. Dissolve 0.7456g potassium chloride into 1 litre pure water to be 0.01N potassium chloride standard solution (1,413µS/cm).

5.2 Calibration mode:

Parameter	Calibration method	Preset temperature compensation method
	Set CELL CONSTANT	(Lin) 2.00% adjustable
Conductivity	Std 1P Prepare 5 selectable kinds of Buffer	According to standard solution default
	Std 5P Prepare 5 selectable kinds of Buffer	According to standard solution default
Resistivity	Set CELL CONSTANT	(nLin) non-liner temp. compensation

5.3 Entry of calibration mode

Temperature coefficient produces a great influence on measurement of conductivity, so it is suggested to control the temperature of standard solution within $25 \pm 3^{\circ}$ C during the calibration.

5.3.1 Cell constant calibration:

1. When the calibration parameter is set as cell constant calibration mode (see ch4.4), press **CAL** key to enter the electrode coefficient set-up page, and then press **ENTER** key to enter the next page.



- Press upward key or downward key to select cell constant default 0.475, 0.100 or 0.010, 10.00^{1/cm} (resistivity only includes 0.010, 0.050, .0100^{1/cm}), in order to select a suitable cell constant value. After selecting the suitable default, press ENTER key to enter the next page.
 - 3. Then, the cell constant value begins to twinkle. Press upward key or downward key to adjust the cell constant value until the measurement value in the display equal to the buffer's standard. Then, press **ENTER** key to back to the measurement mode.

5.3.2 Standard solution single-point calibration:

- 1. When the calibration parameter is set as single-point calibration (see ch4.4), press **CAL** key to enter Standard solution single-point calibration, and then press **ENTER** key to enter the calibration solution selection page.
- Show the first known buffer's standard and temperature. Press UP or DOWN key to adjust the known buffer's standard, including five default values 10.00µS/cm, 84.0µS/cm, 1413µS/cm, 12.88mS/cm and 100.0mS/cm. After selecting the suitable default, press ENTER key to enter the next page.
- 3. The symbol Auto Read begins to twinkle, indicating that the instrument begins auto calibration.











4. After the calibration, the display shows the calibrated electrode coefficient automatically. Press **ENTER** key to back to the measurement mode.

5.3.3 Standard solution five-point calibration:

- 1. When the calibration parameter is set as **five-point calibration** (see ch4.4), press **CAL** key to enter **Standard solution five-point calibration**, and then press **ENTER** key to enter the first CAL1 page. Press **MODE** key at any time will allow you intermit the calibration and back to the measurement mode.
- Press upward key or downward key to select calibration point CAL1~5. Their default values are respectively 10.00uS, 84.0uS, 1413uS, 12.88mS and 100.0mS. You can select any calibration point to make single-point, two-point, three-point, four-point, and five-point calibration. After that, clean and dry the sensor before putting it into the buffer solution. Press ENTER key to enter the calibration page.
- 3. The symbol Auto Read begins to twinkle, indicating that the instrument begins auto calibration. After the calibration, the display shows the calibrated electrode coefficient automatically.
- 4. After the calibration, the display shows the calibrated electrode coefficient automatically.
- 5. Press **ENTER** key to enter the second calibration (CAL2) and finish the first point calibration. Press **MODE** key to store the calibration data and the system will exit calibration mode and back to measurement mode.
- 6. Other CAL2~CAL5 calibration applies the same steps as CAL1.
- NOTE: The design of this instrument has no limit of the calibration sequence, so after finishing any point calibration, you can press **ENTER** key to enter the next point calibration, and press **MODE** key to store the calibration data. Then, exit the calibration mode and back to the measurement mode to finish the calibration of any one or a few points.

5.3.4 Error messages

In case of any error messages when the calibration fails, please refer to chapter 9, Trouble Shooting, to solve the problem and restart the calibration.



6. Instruction set for RS-232

The communication between the instrument and the Windows Hyper Terminal is set as follows:

1. Click [Start] on the lower left, and select [All Programs] \rightarrow [Accessories] \rightarrow

Communication \rightarrow **Hyper Terminal**, and the display will show a window of "Hyper Terminal".

- 2. Set icon and online name for Hyper Terminal before entering the "ONLINE" window.
- 3. Select communication terminal (for example, select 「connect COM1」) in 「Use online (N):」 in "ONLINE" window, and then enter the window "COM1 content".
- In 「transmission bit per second (B): 」 in "COM1 content", select baud rate=9600 (It should be consistent with conductivity meter settings), data bit (D)=8(compulsory), parity check(P)=non-parity, stop bit (S)= 2, flow control (E) =none.
- 5. Click [disconnet]
- 6. Click [file] \rightarrow [content], and enter the window "XXX content".
- 7. Click [setting value] \rightarrow [ASCII setting] to enter "ASCII setting" window.
- 8. In "ASCII setting" window, select 「newline at the end of each line (S)」,「Respond to the input characters (E)」,「Add LF at the end of each input line (A)」 and 「newline if exceeding the width of terminal (W)」, and then leave "ASCII setting" window by pushing 【confirm】.
- 9. Push [confirm] to leave "XXX content" window.
- 10. Push [call] to link the communication.

Note:

COM 1 content setup				
ASCII Setting				
Baud rate	9600			
Data bit	8 (compulsory)			
Parity	None			
Stop bit	2			
Flow control	None			

7. Modbus protocol and instructions

7.1 Introduction

The SC-110 meter applies standard MODBUS protocol, and supports RTU transmission mode, even parity, and it allows Modbus protocol compatible Master or a controlling software from the third party to proceed in data transmission. With the ability of the meter, you only need to add a set of central communication controlling display software to establish a set of controlling communication network system. It will be convenient to the system monitoring and controlling and for data collection.

The default setting:

ID: 01 Transmission frame mode: ASCII Parity type: Non-parity Baud rate: 9600

7.2 Modbus rule

- 1. All the RS-232 communication loops follow Master/Slave way. According to this way, data can be transmitted from a Master (ex: PC) to a Slave (ex: SC-110).
- 2. The master can initialize and control all the information transmission within the RS-232 loop.
- 3. All the communication cannot start from a Slave.
- 4. All the communication within the RS-232 loop is necessary to match the Modbus information frame.
- 5. If the Master or the slave receives the information frame which is included unknown command, the master or the slave does not respond. Which message format is a byte (data frame) composed of a string.

7.3 Data frame format and data command

The communication transmission is in an asynchronous way, and the unit of it is byte (data frame). Each data frame is in an 11 bits (MODBUS RTU) or 10 bits (MODBUS ASCII) sequence data procedure.

BIT	Modbus RTU	Modbus ASCII
Start bit	1 start bit	1 start bit
Data bits	8 data bits	7 data bits
Parity bit	1 bit for even/odd parity, no bit for no parity	
Stop bits	1 stop bit if even/odd parity is used, 2 bits if no parity	

Data frame format (byte):

7.4 SC-110 communication protocol

When a communication command which is sent by a Master is received by the meter (slave), only the corresponding address code of the device accepts the communication command. If not wrong, then the meter performs the corresponding task; then the meter sends the implementation result back to the sender (Master). The returned information includes an address code, a function code which performs actions, data after performing actions, a check code (CRC or LRC) which checks errors. If there is an error, the meter will send corresponding error message.

7.4.1 Message Framing

RTU	

Start	Address	Function	Data	CRC	End
≥3.5 byte time	Address code	Function code	Data field	Check code	≥3.5 byte
	1 byte	1 byte	N bytes	2 bytes	time

Under RTU mode, the max. length of the info data is 256 bytes.

ASCII

Start	Address	Function	Data	LRC	End
:	Address code	Function code 2 bits	Data field	Check code	CRLF
1 bit	2 bits		N bits	2 bits	2 bits

Under ASCII mode, the maximum length of the info data is 513 bits (char).

7.4.2 Address code

Each slave has its one and only address code. That the master sends the address indicates only the slave with the specific address code can receive the message and responds to the master. The slave returns the address to prove where the slave belongs.

Note: The fault setting of address code for SC-110 is 01 and cannot be revised.

7.4.3 Function code

The range of function code is $1\sim255$; the function codes within $1\sim127$ are normal, but the function codes $128\sim255$ are defined as abnormal responding function code. SC-110 only supports some of the function codes. When starting communication, the Master sends the request through function code to tell the Slave to perform which actions. The Slave responds function which is the same as that sent by the Master to indicate the Slave has already responded to the Master and has performed the action.

The function codes which the slaves normally respond are the same as those from the Master. When the slaves respond abnormal message, the highest bit (MSB) of the function codes are set to be 1. It means that there occurs an error when executing the function codes.

Function code	Definition	Operation		
01H	Reads the ON/OFF status of discrete	Read the status of one or more discrete		
Read Coil Status	cons in the slave	contacts		
03H	Read the binary contents of holding	Read one or more data register value		
Read Holding Register	registers in the slave			
05H	Writes a single coil to either ON or	Set the status of discrete individual		
Write Single Coil	OFF.	contacts to the specified address		
06H	Writes a value into a single holding	Set a single 16 bit data register		
Write Single Register	register	Set a single 10-bit data legister		
0FH	Writes each coil in a sequence of	Set the status of multiple discrete		
	agila to gither ON or OFF	continuous contacts to the specified		
Write Multiple Coils	cons to entier ON of OFF	address		

The following table lists the function codes supported by SC-110

10H Write Multiple Registe	Writes values into a sequence of holding registers	Set multiple sequence of 16-bit data register
08H Diagnostics	Diagnostics function	The assessment of network communication capability

7.4.3.1 Function code 01H (Read Coil Status)

The function code is for reading the consecutive contact state discrete volume from the slave. 01H does not support broadcast mode.

Sending format:

Ty	pe	RTU	ASCII	Example (RTU)	
Function		1 byte	2 Characters	01H	Read the contact state discrete volume
DATA	The starting address of the contact number	2 bytes	4 Characters	0070H	The contact numbers to read the starting address is 0070H
DATA Field	The number of contact	2 bytes	4 Characters	0003H	Start to read three consecutive contact state discrete volume from 0070H

Normal response format:

Ту	pe	RTU	ASCII	Example (RTU)	
Function		1 byte	2 Characters	01H Res cod	ponse to the function e
	Number of bytes	1 byte	2 Characters	01H The byte	e number of data es
DATA Field	The contact state discrete value	N bytes	2*N Characters	03H Set disc If th is n will the	of state of contact crete value. Ex: 03H. ne number of contacts ot a multiple of 8, it l be times 10(add 0 in back) until it is of 8.

For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding address list of function code 01H.

7.4.3.2 Function code 03H (Read Holding Registers)

The function code is for reading the consecutive 16 bits register data of the slave. 03H does not support broadcast mode.

Ту	pe	RTU	ASCII		Example (RTU)
Function		1 byte	2 Characters	03H	Read register data
DATA	The starting address of the register	2 bytes	4 Characters	0004H	Read the starting address of register address 0004H
Field	Number of registers	2 bytes	4 Characters	0003H	Start from 0004H to read three consecutive 16 bits register data

Sending format:

Normal response format:

Туре		RTU	ASCII		Example	(RTU)
Function		1 byte	2 Characters	03H Response Function C		Function Code
DATA	Number of bytes	1 byte	2 Characters	06H	The numb	per of data bytes
Field	Register value	N bytes	2*N Characters	0306050	40303Н 16-b	Returns three oit register data

For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding address list of function code 03H.

Note: Due to all the floating-point data of the SC-110 are in a 32 bits IEEE 754 format. The high bit is in front, and the low bit in post.

7.4.3.3 Function code 05H (Write Single Coil)

The function code is for setting single state discrete contacts to the slave. 05H does not support broadcast mode.

Sending format:

	Туре	RTU	ASCII	Example (RTU)
Function code		1 byte	2 Characters	05H set the status of discrete individual contacts
	Contact number address		4 Characters	0076H set single contact number address as 0076H
Field	The discrete status of contact	2 bytes	4 Characters	0000H or FF00H: 0000H indicates OFF, FF00H indicates ON

Normal response format is the copy of the sending format. It returns when successfully set the discrete status of contact. For abnormal response please refer to abnormal data format.

7.4.3.4 Function code 06H (Write Single Register)

The function code is for setting single register value to the slave. 06H supports broadcast mode.

Sending format:

Туре		RTU	ASCII	Example (RTU)	
Function code		1 byte	2 Characters	06H	Set the value of a single register
DATA Field	Register address	2 bytes	4 Characters	000AH	Set a single register address as 000AH
	Register data	2 bytes	4 Characters	0003H	Set a single register address as 0003H

Normal response format is the copy of the sending format. It returns when successfully set the discrete status of contact. For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding address list of function code 06H.

7.4.3.5 Function code 0FH(Write Multiple Coils)

The function code is for setting continuous state discrete contacts to the slave. 0FH does not support broadcast mode.

Туре		RTU	ASCII	Example (RTU)	
Function code		1 byte	2 Characters	0FH	Set the status of discrete individual contacts
DATA Field	Contact starting address	2 bytes	4 Characters	0076H	Set multiple consecutive contacts starting address as 0076H
	The number of contact	2 bytes	4 Characters	0003H	To indicate to set three numbers of contact
	The number of bytes	1 byte	2 Characters	01H	To set the number of information byte of contact discrete status
	The discrete status of multi-contact	N bytes	2*N Characters	03H	Set of discrete contact state information. If the number of contacts is not a multiple of 8, it will be times 10(add 0 in the back) until it is of 8

Sending format:

Normal response format is the copy of the sending format. It returns when successfully set the discrete status of contact. For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding address list of function code 0FH.

7.4.3.6 Function code 10H (Write Multiple Registers)

The function code is for setting multi-consecutive register value to the slave. 10H supports broadcast mode.

Туре		RTU	ASCII	Example (RTU)	
Function		1 byte	2 Characters	10H To set multi-consecutive register data	
	Register starting address	^{er} g 2 bytes 4 Characters		0007H Set multiple consecutive contacts starting address as 0007H	
DATA Field	The number of register	2 bytes	4 Characters	0003H To indicate to set three consecutive 16 bits register numbers	
	The number of bytes	1 byte	2 Characters	06H To set the information bytes number of multi-consecutive register	
	Register data values	N bytes	2*N Characters	030505030303H To set the input of three 16 bits data 0305H, 0503H, 0303H	

Sending format:

Normal response format is the copy of the sending format. It returns when successfully set the discrete status of contact. For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding address list of function code 10H.

7.4.3.7 Function code 08H (Diagnostics)

The function code, 08H, is for diagnostic function. It can be counted packets of every state to evaluate transmission capacity of RS-232 communication.

The function code 08H provides a series of sub-function code. The SC-110 supports sub-function code 0A-12H. The function code 08H does not support broadcast mode.

Ту	уре	RTU	ASCII	Example (RTU)	
Function code (CS)		1 byte	2 Characters	08H Diagnostic function	
DATA	Sub- function code	2 bytes	4 Characters	000AH	Clear counters
Field	Data	2 bytes	4 Characters	0000H	Sub-function code 0A-12H is fixed to 0

Sending format:

T	ype	RTU	ASCII	Example (RTU)	
Func	tion code	1 byte	2 Characters	08H	Response function code
	2 bytes	2 bytes	4 Characters	000AH	Response sub-function code
DATA Field	2 bytes	2 bytes	4 Characters	0000H	Back to the counter value Only sub-function code 0A is able to copy the data and to send the information

Normal response format:

For abnormal response please refer to abnormal data format.

Important: To see the detailed operation address of SC-110, please refer to the corresponding diagnostics function of function code 08H.

7.4.4 Data field

Data field varies with the function code. Whether address or register, the information is high byte first and low byte in the post. The byte sets of the ASCII are twice than those of RTU mode.

7.4.5 Check code

Check code is used to detect whether data frame is error or not. If the data frame is error, the data does not work. It ensures the safety and efficiency of the system. RTU mode uses CRC (loop redundant) to check. ASCII mode uses (Longitudinal Redundancy) LRC to check.

7.4.6 Abnormal procedure

SC-110 will response abnormal information frame when it detects error which except verification error and length of byte error. The minimum byte(LSB) of function code is the same as it of the master, however, the maximum byte (MSB) of function code is 1. It means that the function code which is responded by remote equipment occur has errors. is based on the function code which is sent by master add 128.

Abnormal response frame:

Function code	Abnormal code
MSB: 1	01 or 02 or 03 or 04

Abnormal code 01: illegal function code

The SC-110 does not support the function code received.

Abnormal code 02: illegal data address

The SC-110 do not support the designated data address.

Abnormal code 03: illegal statistics value

The data which is input to designated address of SC-110 is illegal value.

Abnormal code 04: abnormal data input

Failed to input data to SC-110, and it result to unrecoverable error.

7.5 Communication connection

The RS-232 communication port of the meter features with electronic isolation protection, lightning protection, and to provide internal independent ground solution. It is allowed to use normal twisted-pair (segregation double-stranded twisted pair cable) cable connections. All devices are in contact with a double-stranded, and then all together, and another line will be connected with all the negative contacts, and the isolated shield wire must be connected to GND. When we talk about communication in the laboratory, the stand-alone master-slave communication is relatively simple. Hence, it is allowed to consider using the normal cable instead. However, there should be strictly in accordance with the requirements of industrial engineering construction.

7.6 MODBUS name and address table

Logic	Item	Number	Information	Description of	Default value	Note
address	Item	of Byte	type	data transmission	Delaunt value	NOIC
0001H	Equipment's ID	2	USHORT	1-247	1	
0002H	Meter model	6	USHORT	ASCII Code	SC110	
000511	Communication	2	UCHODT	0: RTU	0	
0005H	protocol	2	USHORT	1: ASCII	0	
				0:2400		
00061	Serial transmission	2	USHOPT	1:4800	2	
00001	speed (Baud rate)	2	USHOKI	2:9600	5	
				3: 19200		
				0: None		
0007H	Parity	2	USHORT	1: Even	1	
				2: Odd		
0008H			USHORT	Second		
0009H			USHORT	Minute		
000AH	Real-time clock* 12	USHORT	Hour	2000-01-01,		
000BH		USHORT	Day	00:00:00		
000CH			USHORT	Month		
000DH			USHORT	Year		
000EH	Factory reserved					
				0: MTC		
000FH	Temperature mode*	2	USHORT	1: PTC	0	
	•			2: NTC		
0010H						
0011H						
0012H	Factory reserved					
0014H	-					
0015H						

Function Code: 03H, 06, 10H Modbus response (setup parameter)

Function code: 03H Modbus response (measurement parameter)

Logic address	Item	Number of Byte	Information type	Description of data transmission	Default value	Note
0031H	Number of measurement channels	2	USHORT	SC-110 only has one channel	1	
0032H	Sign byte	6	CHAR	μS/cm		ASCII code
0035H	Residual chlorine /Ozone measurement	4	FLOAT	Cond/Res/Salt /TDS measurement		Data affected by sign byte

0037H	Temperature measurement	4	FLOAT	Temperature measurement	
0039H- 0050H			Factory	reserved	

Function code: 08H Correlated diagnostic function

Sub function code	Name of Counter	Note
0AH	Clear all the counters	Clear Counters and Diagnostic Register
0BH	Return Bus Message Count	The response data field returns the quantity of message that the slave has detected on the communications system since its last restart or last clear counters operation, or being powered-up. It counts whether the address comes from the remote equipment or not,
0CH	Return Bus Communication Error Count	The response data field returns the quantity of CRC errors encountered by the slave since its last restart, clear counters operation, or power-up. It counts whether the address comes from the remote equipment or not,
0DH	Return Bus Exception Error Count	The response data field returns the quantity of Modbus exception responses returned by the slave since its last restart, clear counters operation, or power-up. It only counts when the address comes from the remote equipment.
0EH	Return Slave Message Count	The response data field returns the quantity of messages addressed to the slave, or broadcast, that the slave has processed since its last restart, clear counters operation, or power-up. It only counts when the address comes from that remote equipment.
OFH	Return Slave No Response Count	The response data field returns the quantity of message addressed to the slave for which it returned no response (neither a normal responses nor an exception response), since its last restart, clear counters operation, or power-up. It is said that the counter will calculate the quantity of none-error broadcast messages.
10H	Return Slave NAK Count	The response data field returns the quantity of message addressed to the slave for which it returned a Negative Acknowledge (NAK) exception response since its last restart, clear counters operation, or power-up.
11H	Return Slave Busy Count	The response data field returns the quantity of message addressed to the slave for it returned a Slave Device Busy exception response, since its last restart, clear counters operation, or power-up.
12H	Return Bus Character overrun Count	The response data field returns the quantity of messages addressed to the slave that it could not handle due to a character overrun condition, since its last restart, clear counters operation, or power-up.

8. Optional accessories

Product Name	Order Number
4-electrode conductivity probe, $C = 0.475$	8-243
Ultra pure water resistivity probe, $C = 0.01$	8-221-01
Pure water resistivity probe, $C = 0.05$	8-222-01
3/4"NPT bypass chamber for 8-22X-01 resistivity probe	8-TF-02
Cable for 8-22X-01 resistivity probe, 1.15M	8-101D
Cable for 8-22X-01 resistivity probe, 3M	8-103D
RS-232 Connection Cable for PC	8-30
Carrying Case for SC-110	8-32-1
Conductivity Calibration Solutions, 1,413µS/cm, 500ml	8-EC1413-2
Conductivity Calibration Solutions, 84µS/cm, 500ml	8-EC84-2

9. Trouble shooting

Failure phenomenon	Possible reason	Disposition
ERR1	The electrode coefficient bears a	Replace new standard solution for another
	large deviation when use the	calibration. For others, please refer to the
	standard solution in calibration.	disposition of measurement deviation in
		the following.
ERR2	The measurement is unstable when	Please refer to the disposition of unstable
	use the standard solution in	measurement in the following.
	calibration.	
Measurement deviation	The electrode is polluted.	Clean and wash the electrode
	The electrode coefficient (C) is	Modify the electrode coefficient
	wrongly set.	
	The temperature coefficient (TC) is	Modify the temperature compensation
	wrongly set.	coefficient.
	The electrode is damaged or	Readjust the coefficient
	broken-down.	
	There are air bubbles on the	Churn up to remove the bubbles.
	electrode.	
	The battery is low.	Replace with new battery.
	Instrument fault	Send to the original factory for repair.
Unstable measurement	The electrode is disturbed by the	Churn up to remove the bubbles.
	bubbles.	
	The electrode is polluted.	Clean and wash the electrode.
	Electrode fault	Replace the electrode.
	The battery is low.	Replace with new battery.
	Instrument fault	Send to the original factory for repair.
The conductivity value is	The electrode plugs is not inserted	Check the electrode plugs and reinsert it
zero or the special	into the required place.	into the required place.
resistivity value is out of	The electrode is broken.	Replace the electrode.
the range.		
Temperature deviation	Not soak deeply enough.	Deepen the soaked length to ensure the
		accuracy of temperature.
	Deviation of standard comparison	Use standard thermometer comparison to
		correct the difference in temperature
	The temperature sensor is broken	Replace the electrode.
	The battery is low.	Replace with new battery.
	Instrument fault	Send to the original factory for repair.
Temperature is wrongly	The electrode plugs is not inserted	Check the electrode plugs and reinsert it
displayed.	into the required place.	into the required place.
	The temperature sensor is broken	Replace the electrode.
	Instrument fault	Send to the original factory for repair.

10. Maintenance

10.1 Maintenance of instrument machine

Please keep it in dry and ventilated environment. Use wrong soft wet cloth to wipe the generally polluted surface. For pollution by grease, please remove the stain with waxed soft wet cloth. Solvent is strictly prohibited.

10.2 Battery replacement

- a. When the display shows Lo Bat during the instrument works, please replace with alkali dry cells.
- b. First power off the instrument, and then use cross head screwdriver to dismantle the four screws on the back of the instrument.
- c. Open the back cover and remove the battery protection cover.
- d. Take out the old battery and put in 4 new AA batteries in correct directions.
- e. Reinstall the battery cover and back cover, and lock the 4 screws to finish the battery replacement.
- f. Caution! The battery replacement should be finished within 20 seconds; otherwise the date data will disappear and should be reset.

10.3 Maintenance of sensor

The substances on the measuring side of sensor will influence the measurement value and cause deviation. Please clean and maintain it regularly. Use soft brush or cotton stick to clean the graphite surface of sensor.

Pollution types	Cleaning methods	
Pollution by grease	Use mild degreasing agent to clean the electrode.	
Pollution by dirt	Soak in 2~5% HCl for about 5 minutes.	



SUNTEX INSTRUMENTS CO., LTD.

13F, No. 31, Lane 169, Kangning St., Xizhi Dist., New Taipei City, Taiwan (R.O.C.) Tel: 886-2-2695-9688 Fax: 886-2-2695-9693 e-mail: suntex@ms1.hinet.net

www.suntex.com.tw/en

5040E SC-110/ Technical data subject to alternations/ Quality Systems ISO 9001/201612