




ThermoFisher
S C I E N T I F I C


X선 형광분석 소개



Introduction to **X-Ray** Fluorescence

How the XRF Works



- Each individual element produces its own set of characteristic x-rays; the basis for **qualitative** analysis
- By counting the number of characteristic x-rays of a given element we can determine its concentration; the basis for **quantitative** analysis






엑스선 입사
(광원)

엑스선 형광
(샘플)

각 원소별 분석
(검출기)

정성 및 정량
분석
(내장 PC)


Technology Basics



X-ray Energy Reference



H Hydrogen 1																	He Helium 2						
Li Lithium 3	Be Beryllium 4																	B Boron 5	C Carbon 6	N Nitrogen 7	O Oxygen 8	F Fluorine 9	Ne Neon 10
Na Sodium 11	Mg Magnesium 12																	Al Aluminum 13	Si Silicon 14	P Phosphorus 15	S Sulfur 16	Cl Chlorine 17	Ar Argon 18
K Potassium 19	Ca Calcium 20	Sc Scandium 21	Ti Titanium 22	V Vanadium 23	Cr Chromium 24	Mn Manganese 25	Fe Iron 26	Co Cobalt 27	Ni Nickel 28	Cu Copper 29	Zn Zinc 30	Ga Gallium 31	Ge Germanium 32	As Arsenic 33	Se Selenium 34	Br Bromine 35	Kr Krypton 36						
Rb Rubidium 37	Sr Strontium 38	Y Yttrium 39	Zr Zirconium 40	Nb Niobium 41	Mo Molybdenum 42	Tc Technetium 43	Ru Ruthenium 44	Rh Rhodium 45	Pd Palladium 46	Ag Silver 47	Cd Cadmium 48	In Indium 49	Sn Tin 50	Sb Antimony 51	Te Tellurium 52	I Iodine 53	Xe Xenon 54						
Cs Cesium 55	Ba Barium 56		Hf Hafnium 72	Ta Tantalum 73	W Tungsten 74	Re Rhenium 75	Os Osmium 76	Ir Iridium 77	Pt Platinum 78	Au Gold 79	Hg Mercury 80	Tl Thallium 81	Pb Lead 82	Bi Bismuth 83	Po Polonium 84	At Astatine 85	Rn Radon 86						
Fr Francium 87	Ra Radium 88																						
			La Lanthanum 57	Ce Cerium 58	Pr Praseodymium 59	Nd Neodymium 60	Pm Promethium 61	Sm Samarium 62	Eu Europium 63	Gd Gadolinium 64	Tb Terbium 65	Dy Dysprosium 66	Ho Holmium 67	Er Erbium 68	Tm Thulium 69	Yb Ytterbium 70	Lu Lutetium 71						
			Ac Actinium 89	Th Thorium 90	Pa Protactinium 91	U Uranium 92	Np Neptunium 93	Pu Plutonium 94	Am Americium 95	Cm Curium 96	Bk Berkelium 97	Cf Californium 98	Es Einsteinium 99	Fm Fermium 100	Md Mendelevium 101	No Nobelium 102	Lr Lawrencium 103						




Key to Energy values



Light Elements

57-71
89-103

XRF Spectrum



← General Metals

#535 5.9 sec

SS-316 0.25
Excellent

Ele	%	±2σ
Fe	68.999	0.446
Cr	16.071	0.239
Ni	10.114	0.312
Mo	2.085	0.040
Mn	1.146	0.182
LEC	0.500	0.000
Cu	0.374	0.087


Application for HHXRF




Energy Markets /
Positive Material
Identification (PMI)



Fabrication and
Metals QC/QA



Scrap Recycling



Weld analysis



ThermoFisher
SCIENTIFIC

514659550

XL5 Coatings Mode

Mathieu Bauer, Senior Application Scientist
SEATW Sales Meeting, Bali May 9-11 2018



Coatings Mode Overview



- Designed to determine the coat weight or the coating thickness of metal coatings
- Enables testing of large or irregularly shaped samples, as well as small-diameter wiring or tubing and direct testing of coated samples at the plating line
- Enables analysis of up to four coating layers
- Layer can be pure element, alloy or compound
- Substrate can be pure element, alloy, plastic or wood.
- Optional 3mm small spot analysis
- Element suite - 30 elements including:
 - Pb, **Hg**, Au, Pt, **W**, Sb, Sn, Cd, Ag, Pd, **Rh**, **Ru**, Mo, Nb, Zr, **Y**, Zn, Cu, Ni, Co, Fe, Mn, Cr, V, Ti, S, P, Si, Al, Mg
- Multi point standardization feature enables user to easily optimize calibration with known standards.



Markets and applications

Coatings generally applied for decorative purpose or for improving corrosion wear and heat resistance solderability or electrical conductivity.

- Metal Finishing (40% of total avail. market)
 - Automotive (Cr/Ni/Cu/Plastics, ZnNi/LAS, passivation coatings)
 - Aerospace (ZnNi/LAS, Cd/LAS, etc.)
 - Fasteners (NiP/LAS, Zn/LAS, etc.)
 - Fixtures (Cr/Ni/Brass, Cr/Ni/Cu/ABS etc.)
 - Cutting Tools (WC/TS etc.)
 - Electrical appliances (Zn/Fe NiP/LAS etc.)
 - Accessories, apparel (Au/Ni/Cu/Zn etc.)
 - Jewelry & Watchmaking (Rh/Brass etc.)
 - Electrification products (Ag/Cu etc.)
- Electronics (60% of total avail. market)
 - PCB
 - Connectors
 - Hard drives



Handheld XRF suitable for metal finishing but not in electronics



Coating Mode Basics

- Using General Metal Mode assumes to analyze homogeneous metals and alloys. Units of the measure is wt.%, ppm, etc.
- Coated materials are heterogeneous: composition at the surface is different from composition of the core. General Metal Mode not accurate

Coatings Mode conceived to measure coat weight or layer thicknesses applied on a substrate, not to measure the composition of layer or substrate

- Units designating quantity of coating: applicable to surface, not to bulk composition.
 - Layer Thickness: unit of length e.g. μm
 - Coat Weight: mass per unit of area, e.g. g/m^2

In coatings analysis, the user needs to define substrate material and the sequence of layers



Example



← General Metals

#488 3.4 sec

No Match

Ele	Min	Meas	Max
Ni		64.25	
Cu		29.79	
Zn		3.55	
Cr		1.81	
Pb		0.48	
Fe		0.06	
Below LOD			4σ

← Cr/Ni/Brass




#489 7.3 sec

Layer	μm	±2σ
2: Cr	0.15	0.00
1: Ni	10.38	0.06
CuZn33Pb		Substrate

Coatings Mode Profile. Get Started!

- Unlike General Metals and other Modes, Coatings Mode can't operate without defining the substrate and the sequence of layers


Create a profile to describe coated material




Create a profile

Example of Cr/Ni/Brass substrate

- Assign name to the profile
- Select number of layers




Enter Profile Name



Select the number of layers (up to 4)

Define sequence and nature of layers


- Layer 1: the closest to the substrate, layer with highest # (up to 4) at the surface
- Select layer material
- Density of the layer (default value of pure metal or alloy) can be edited in case an element is mixed with some compound
- Slope of the layer can be adjusted to match values of available reference samples



Select layer to edit




Layer Edit






Selection of layer material

Layers: pure elements. Ni, Cr




Select among list




Define the material of the substrate


Select type of Material: Alloy




Define Alloy Composition or select library



Nominal or measured chemistry



Pick alloy & avg. chemistry from library



Review of the profile, analysis of sample

Edit Cr/Ni/Brass

Mode
Coatings

Profile Name
Cr/Ni/Brass

Number of Layers
2

Layer 2
Element: Cr
Density: 7.180 g/cm³

Layer 1
Element: Ni
Density: 8.876 g/cm³

Substrate
Alloy: C464NavBs

Display Units
µm

Scan Times & Filters
Max Scan Time: 300
Main: 30

Use Data Field Set
Default Data Fields

Enable Small Spot

Beep Times (1st, 2nd, 3rd)
10, 20, 30 seconds

Save Images with Readings
None

Standardization Data
Not Standardized

Standardize On Next Reading

Pseudo Layers
None

Cancel Save

Results of analysis 2 µm Cr foil/ 10µmNi Foil/ C464 substrate

← Cr/Ni/Brass

#364 10.7 sec

Layer	µm	±2σ
2: Cr	2.21	0.01
1: Ni	9.19	0.05
C464NavBs	Substrate	




Result "Out of the Box"

Measurement units

- Layer Thickness: Ang, μm , mm, μi , mils, thou
- Coat Weight: $\mu\text{g}/\text{cm}^2$ g/m^2 , mg/m^2 , g/cm^2 , mg/cm^2 , g/ft^2 , mg/ft^2 , $\mu\text{g}/\text{ft}^2$, oz/ft^2
- Coat weight mostly used when coating material is not pure metal but compound (metal passivation, etc.)

$$Th = \frac{CW}{d}$$

Th: layer thickness
CW: coat weight
d: density




Selectable units

Capabilities of coatings Mode

Analyzing up to 4 layers

2 μ mCr/10 μ mNi/1.5 μ mCu/
6.3 μ mZn/Fe substrate




← test

#373 9.5 sec

Layer	μ m	$\pm 2\sigma$
4: Cr	2.35	0.01
3: Ni	11.06	0.08
2: Cu	1.49	0.03
1: Zn	6.12	0.16
Fe	Substrate	

Analyzing multiple layers over Plastics

Common decorative coating




← Cr/Ni/Cu/plastic

#535 16.8 sec

Layer	μ m	$\pm 2\sigma$
3: Cr	0.15	0.00
2: Ni	17.57	0.11
1: Cu	17.01	0.35
Plastic	Substrate	

Analyzing thin layers

Passivation coating
38mg/m² Zr/HDG



← Zr HDG

#564 4.3 sec

Layer	mg/m ²	$\pm 2\sigma$
2: Zr	40.24	5.26
1: Zn	82723.41	478.14
Fe	Substrate	

Accurate Determination of Multilayer Coat Weight & Thickness

Capabilities of coatings Mode. Alloy substrates

Determination of Ni Monolayers over Stainless steel

10µm Ni/ SS 304

Layer	µm	±2σ
1: Ni	9.85	0.05
SS-304	Substrate	

SS-304: Cr:17.9%, Ni:8.36%,
Mn:1.6%, Mo:0.29%, Fe:Bal

10µm Ni/SS 316

Layer	µm	±2σ
1: Ni	9.62	0.05
SS-316	Substrate	

SS-316: Cr:16.7%, Ni:10.2%,
Mn:1.6%, Mo:2.14%, Fe:Bal

10µm Ni/SS 310

Layer	µm	±2σ
1: Ni	9.60	0.06
SS-310	Substrate	

SS-310: Cr:24.9%, Ni:19.7%,
Mn:1.6%, Mo:0.39%, Fe:Bal

**Accurate determination of layer thickness over alloy substrates
containing the same element**

Capabilities of coatings Mode. Alloy as a layer

Determination of Alloy layers over pure elements and alloys

10.14µm NiZn/ Fe

Layer	µm	±2σ
1: ZnNi6	10.20	0.04
Fe	Substrate	

NiZn: Ni:6% Zn:94%

8.0 µm NiP/Cu

Layer	µm	±2σ
1: NiP	7.84	0.04
Cu	Substrate	

Electroless Ni plating
NiP: 88%Ni 12%P

20 µm NiP/Kovar


Layer	µm	±2σ
1: NiP	18.58	0.11
Kovar	Substrate	

Electroless Ni plating
NiP: 92%Ni 8%P
Kovar: 29%Ni%, 17%Co, Bal Fe

Accurate determination of thickness of alloy layers over pure metals or alloys substrates

Standardization

- Check Standardization on next reading, enter # of standards (max.5)
- When prompted, enter reference value for each standard
- Perform the measurement
- After last measurement, linear correction is calculated and applied



Kovar: Ni:29%, Co:17%, Fe:Bal

Edit Ni Kovar

Mode
Coatings

Standardization Data
13:10 03/14/2017

Standardize On Next Reading

Readings to Standardize
4

Pseudo Layers
None

Cancel Save

Enter 1 of 4

Layer 1: Ni
10.240

Cancel OK

ni kov

#240 15.8 sec

Layer	μm	$\pm 2\sigma$
1: Ni	9.53	0.04

Kovar Substrate

Cancel OK

Standardization Data

Layer 1
Slope: 1.025, Intercept: 0.272, R^2 : 0.999865

Cancel OK

Pseudo Layers: compounds as coating

Define Pseudo Layer as function of existing layers

Example P as P2O5/
Galvanized steel

340mg/m² P2O5/Galvanized

Pseudo Layer

P2O5
L2*2.29

1	2	3	+	-
4	5	6	*	/
7	8	9	()
⊗	0	.	Layer	

Cancel OK

← phosphate

#425 17.0 sec

Layer	mg/m ²	±2σ
2: P	163.54	4.10
1: Zn	41609.93	161.32
Fe		Substrate

Pseudo Layers


P2O5	374.51
------	--------

Determination of Compounds Coat Weight

Infinite Thickness – Saturation

50 μ m Ni /Cu substrate

Layer	μ m	$\pm 2\sigma$
1: Ni	Saturated	
Cu	Substrate	




T: Tube
D: Detector

Saturation: no signal increase at increasing thickness

Infinite Thickness – Buried Layer

15µm Zn/20µmNi/2µmCr/AA6063

Layer	µm	±2σ
3: Zn	14.88	0.46
2: Ni	18.77	2.86
1: Cr	Buried	
AA 6063	Substrate	



T: Tube
D: Detector

Buried layer: signal from substrate & deep layer(s) not detected

Infinite Thickness – examples

Calculated values of infinite thicknesses (80% absorption) for different coating/substrate

Coating	Substrate	Infinite thickness of the coating element (μm).	Infinite thickness of the coating for the substrate element (μm)
Ag	Cu	59	6.0
Au	Ni	8.0	2.9
Au	Cu	8.0	3.5
Cd	Fe	73	4.0
Cr	Fe	22	4.0
Cr	Ni	22	6.2
Cu	Al	29	0.3
Cu	ABS	29	N/A
Ni	Al	26	0.4
Ni	Cu	26	31
Ni	Fe	26	17
Pb	Cu	22	7.7
Pd	Ni	51	4.7
Rh	Cu	48	5.8
Sn	Al	89	0.6
Sn	Cu	89	7.5
Ti	Fe	28	8.7
Zn	Fe	39	18
Zr	Al	79	1.3