Colilert^{*}

Predispensed MPN

06-14784-02

IDEXX

For Technical Support, please call:

North/South America: 207 556 4496/1 800 321 0207

Europe: 00800 1234 3399

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: +81 422 71 5921

Australia: 1800 443 399

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

idexx.com/water

Colilert* Test Kit

Introduction

Colilert* simultaneously detects total coliforms and *E. coli* in water. It is based on IDEXX's patented Defined Substrate Technology* (DST*). When total coliforms metabolize Colilert's nutrient-indicator, ONPG, the sample turns yellow. When *E. coli* metabolize Colilert's nutrient-indicator, MUG, the sample also fluoresces. Colilert can simultaneously detect these bacteria at 1 cfu/100 mL within 24 hours even with as many as 2 million heterotrophic bacteria per 100 mL present.

Storage

Store at 2–30°C, away from light.

MPN Test Procedure

1. Select the appropriate number of tubes per sample for your MPN test (5,10, etc.).
2. Aseptically fill each Colilert tube with 10 mL of a well-mixed water sample.
3. Cap the tubes tightly.
4. Mix vigorously to dissolve the reagent by repeated inversion. Some particles may remain undissolved. Dissolution will continue during incubation.
5. Incubate inoculated reagent tubes at $35 \pm 0.5^\circ\text{C}$ for 24 hours. *This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.*
6. Read results according to the Result Interpretation table below. To find the concentration of total coliforms or *E. coli* per 100 mL, compare the number of positive tubes per sample set to the standard MPN (Most Probable Number) probability chart as shown.

Result Interpretation

Appearance	Result
Less yellow than the comparator ¹	Negative for total coliforms and <i>E. coli</i>
Yellow equal to or greater than the comparator	Positive for total coliforms
Yellow and fluorescence equal to or greater than the comparator	Positive for <i>E. coli</i>

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Colilert results are definitive at 24–28 hours. In addition, positives for both total coliforms and *E. coli* observed before 24 hours and negatives observed after 28 hours are also valid.

MPN Index and 95% Confidence Limits for Various Combinations of Positive and Negative Results When Five 10-mL Portions are Used²

No. of Tubes Giving Positive Reaction Out of 5 of 10 mL Each	MPN Index/Lower	95% Confidence Limits (Approximate)	
		Lower	Upper
0	< 2.2	0	6.0
1	2.2	0.1	12.6
2	5.1	0.5	19.2
3	9.2	1.6	29.4
4	16.0	3.3	52.9
5	> 16.0	8.0	Infinite

MPN Index and 95% Confidence Limits for Various Combinations of Positive and Negative Results When Ten 10-mL Portions are Used²

No. of Tubes Giving Positive Reaction Out of 10 of 10 mL Each	MPN Index/100 mL	95% Confidence Limits (Approximate)	
		Lower	Upper
0	< 1.1	0	3.0
1	1.1	0.03	5.9
2	2.2	0.26	8.1
3	3.6	0.69	10.6
4	5.1	1.3	13.4
5	6.9	2.1	16.8
6	9.2	3.1	21.1
7	12.0	4.3	27.1
8	16.1	5.9	36.8
9	23.0	8.1	59.5
10	> 23.0	13.5	Infinite

Procedural Notes

- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures. For example, samples run in other countries are incubated at $36 \pm 2^\circ\text{C}$ for 24–28 hours.
- Colilert can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*² MPN tables should be used to find Most Probable Numbers (MPNs).
- If a water sample has some background color, compare inoculated Colilert sample to a control blank of the same water sample.
- If sample dilutions are made, multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- Colilert is a primary water test. Colilert performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- In samples with excessive chlorine, a blue flash may be seen when adding Colilert. If this is seen, consider sample invalid and discontinue testing.
- Aseptic technique should always be followed when using Colilert. Dispose of in accordance with Good Laboratory Practices.

Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Colilert:
 - A. IDEXX-QC Coliform and *E. coli*³: *Escherichia coli*, *Klebsiella pneumoniae*, and *Pseudomonas aeruginosa*.
 - B. Quanti-Cult**⁴ *Escherichia coli*, *Klebsiella pneumoniae* and *Pseudomonas aeruginosa*.
 - C. Fill three sterile vessels with 10 mL of sterile nonbuffered oxidant-free water and inoculate with a sterile loop of ATCC⁵ strains, *Escherichia coli* ATCC 25922 or 11775, *Klebsiella pneumoniae* ATCC 31488 and *Pseudomonas aeruginosa* ATCC 10145 or 27853.
2. Follow the MPN Test Procedure above.
3. Results should match the Result Interpretation table above.

1. IDEXX MPN Comparator, catalog W102

2. Eaton AD, Clesceri LS, Rice EW. *Standard Methods for the Examination of Water and Wastewater*. Washington, DC: American Public Health Association; 2005.

3. IDEXX-QC Coliform and *E. coli*—IDEXX Catalog #UN3373-WQC-TCEC

4. Quanti-Cult cultures—IDEXX catalog # WKIT-1001

5. American Type Culture Collection 1-800-638-6597 atcc.org

*Colilert, Defined Substrate Technology and DST are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. Quanti-Cult is a trademark or registered trademark of Remel Inc.

Patent information: idexx.com/patents.

© 2013 IDEXX Laboratories, Inc. All rights reserved.

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