

Specifications	Mi151
Range pH	-2.00 to 16.00 pH
mV	±699.9 mV / ±1999 mV
Temp	-20.0 to 120.0°C / -4.0 to 248.0°F
Resolution pH	0.01 pH
mV	0.1 mV / 1 mV
Temp	0.1°C (0.1°F)
Accuracy pH	±0.01 pH
(@ 20°C) mV	±0.2 mV / ±1 mV
Temp	±0.4°C / ±0.8°F
Typical EMC pH	±0.02 pH
Deviation mV	±0.2 mV / ±1 mV
Temp	±0.4°C / ±0.8°F
pH Automatic Calibration	1 or 2 point-calibration, with 7 memorized buffers
Offset Calibration	±1 pH
Slope Calibration	from 80 to 108%
Temperature Compensation	automatic, from -20.0 to 120.0°C / -4.0 to 248.0°F or manual, without temperature probe
pH Electrode	MA917B/1 (included)
Temperature Probe	MA831R (included)
Environment	0 to 50°C / 32 to 122°F; max RH 95%
Input Impedance	10 <sup>12</sup> Ohm
Power supply	12 VDC power adapter (included)
Dimensions	230 x 160 x 95 mm
Weight	0.9 kg

# Accessories

MA9001 pH 1.68 buffer solution, 230 mL bottle MA9004 pH 4.01 buffer solution, 230 mL bottle MA9006 pH 6.86 buffer solution, 230 mL bottle MA9007 pH 7.01 buffer solution, 230 mL bottle MA9009 pH 9.18 buffer solution, 230 mL bottle MA9010 pH 10.01 buffer solution, 230 mL bottle MA9012 Refilling solution for double junction electrode, 230 mL bottle

MA9015 Electrode storage solution, 230 mL Electrode cleaning solution, 230 mL MA9016

MA831R Temperature probe









12 VDC Adapter, 220 V 12 VDC Adapter, 110 V Electrode Holder MA9311 MA9315

MA917B/1 Glass body, double junction refillable pH electrode Platinum ORP electrode with 1 m MA921B/1

cable (will be replaced by SE300) SE300 Platinum ORP electrode with 1 m cable

### Glass Electrode & **Temperature** Probe

Choose from our wide selection of pH and ORP electrodes at pages 6 and 45.

#### Custom dual level LCD

Large and easyto-read Custom dual level LCD Display with simultaneous readings and with user-friendly icons.



## Ordering Information

Mi151 is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
- MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution
- M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution
  M10016 Sachet Electrode Cleaning Solution
- Graduate Pipet
- 12 VDC Adapter
- · Instruction manual